## Legionella Biofilm Test

...delivering microbiology results in minutes instead of days

Pinpoint contamination
Find *Legionella* in their breeding ground
Fast, on-site results

The Legionella Biofilm Test<sup>™</sup>, is a breakthrough technology from hydrosense<sup>®</sup> for the rapid detection of *Legionella* bacteria– delivering test results in minutes instead of days.

The test kit contains all the items required to perform 5 tests – including swabs, recovery buffer, and re-suspension tubes.

- **Pinpoint the contamination** in a system to allow targeted cleaning.
- **Demonstrate cleanliness** of an area e.g. after cleaning cooling tower pack.
- A clean water sample does not equate to a clean system: proliferation of *Legionella* is greater in a biofilm than in water.
- Highly selective can detect *Legionella* even in samples with very high background flora.
- Simple to use and fast.
- Results are **easy to read** and obtained within 25 minutes.
- The test allows **immediate understanding** and response to contaminated systems.
- Use can reduce risk of outbreaks and litigation.
- Regular use will **increase confidence** in your *Legionella* control.



Legionella Biofilm Test™ Product code 100144

hydrosense


The hydrosense<sup>®</sup> Legionella Field Test<sup>™</sup> uses an immunochromatographic assay, to detect the presence of cell surface antigens from *Legionella* bacteria in a sample. The presence of antigen causes the "Test Line" to turn red in colour. A "Control Line" is included which should always turn red if the test has been performed correctly.

Springfield, Linlithgow EH49 7NX, United Kingdom. T:+44 (0)1506 840888. F:+44 (0)1506 840777.

www.hydrosense.biz

TECH0020v1.1

Product details	
Bacteria detected	Legionella pneumophila serogroup 1
Detection method	Lateral flow immunochromatographic assay
Lower limit of detection	Up to 200 CFU per area swabbed
Specificity	Tested against a wide range of waterborne bacteria, including Acinetobacter calcoaceticus, Aeromonas hydrophila subsp. hydrophila, Bacillus subtilis, Burkholderia cepacia, Citrobacter freudii, Citrobacter koseri, Escherichia coli, Enterobacter cloacae, Klebsiella oxytoca, Pseudomonas aeruginosa, Pseudomonas fluorescens, Pseudomonas putida, Pseudomonas stutzeri, Ralstonia pickettii, Raoultella terrigena, Streptococcus pyrogenes, Yersinia ruckeri. No cross reactions observed at concentrations $\leq 10^9$ CFU/L. Staphylococcus aureus has been observed to cross react with the test at concentrations $> 10^8$ CFU/L. Legionella pneumophila serogroups 4 & 7 cross react with the test at concentrations $> 10^9$ CFU/sample. The test does not react with any other L. pneumophila serogroups or L. species tested.
Test volume	100 μL Exact volume disposable pipette provided
Suggested sample area	10 cm <sup>2</sup>
Biocides	The product has been tested with a range of common biocides and biodispersants used in <i>Legionella</i> control, at their normal maximum operating concentration. Active ingredients tested include: glutaraldehyde, didecyl-dimethyl-ammonium chloride, isolthazoles, dibromonitrilopropionamide, ionic and non-ionic surfactants. The test gives a false positive if used with polymeric biguanide and THPS.
Operator time	1-2 minutes
Development time	25 minutes
Shelf life	18 months from manufacture when stored at room temperature (<30 $^{\circ}\mathrm{C}$ ) in original packaging



hydrosense®

Springfield, Linlithgow EH49 7NX, United Kingdom. T:+44 (0)1506 840888. F:+44 (0)1506 840777.

www.hydrosense.biz

TECH0020v1.1