Industrial Legionella Test

...delivering microbiology results in minutes instead of days

- Fast results
- High sensitivity

In-line samplingIdeal for industrial samples

ΤМ

The Industrial Legionella Test[™], is a breakthrough technology from hydrosense[®] for the rapid detection of *Legionella* bacteria in industrial water – delivering test results in minutes instead of days.

The test kit contains all the items required to perform 5 tests – including filters, recovery buffer, collection vessel and sample point.

- These tests are optimised for use with **industrial** water systems, **cooling towers** and whirlpool **spas**.
- This filtration step needs **no pumping**, just connect the sample point to the water system and use the system pressure to drive the filtration.
- Single use hollow fibre filtration modules.
- Pre-filled syringe containing recovery buffer, no need to measure or dilute solutions.
- The test is **simple to use** and **fast**, with no interpretation required.
- Results are **easy to read** and obtained within 30 to 35 minutes.
- The test allows **immediate understanding** and response to the water conditions.
- Use can **reduce risk** of outbreaks and litigation.
- Regular use will **increase confidence** in your *Legionella* control.



Industrial Legionella Test™ Product code 100182

hydrosense


The hydrosense[®] Legionella Field Test[™] uses an immunochromatographic assay, to detect the presence of cell surface antigens from *Legionella* bacteria in a sample. The presence of antigen causes the "Test Line" to turn red in colour. A "Control Line" is included which should always turn red if the test has been performed correctly.

Springfield, Linlithgow EH49 7NX, United Kingdom. T: +44 (0) 1506 840888. F: +44 (0) 1506 840777.

www.hydrosense.biz

TECH0019v1.0

Bacteria detected	Legionella pneumophila serogroup 1
Detection method	Lateral flow immunochromatographic assay
Lower limit of detection	Up to 100 CFU/L (0.1 CFU/mL)
Specificity	Tested against a wide range of waterborne bacteria, including Acinetobacter calcoaceticus, Aeromonas hydrophila subsp. hydrophila, Bacillus subtilis, Burkholderia cepacia, Citrobacter freudil Citrobacter koseri, Escherichia coli, Enterobacter cloacae, Klebsiella oxytoca, Pseudomonas aeruginosa, Pseudomonas fluorescens, Pseudomonas putida, Pseudomonas stutzeri, Ralstonia pickettii, Raoultella terrigena, Streptococcus pyrogenes, Yersinia ruckeri. No cross reactions observed at concentrations $\leq 10^9$ CFU/L. Staphylococcus aureus has been observed to cross react with the test at concentrations > 10 ⁸ CFU/L. Legionella pneumophila serogroups 4 & 7 cross react with the test at concentrations >10 ⁹ CFU/sample. The test does not react with any other L. pneumophila serogroups or L. species tested.
Test volume	100 µL
Suggested sample volum	1e 250 mL
Sample parameters	Sample should filter in less than 10 minutes.
Biocides	The product has been tested with a range of common biocides and biodispersants used in <i>Legionella</i> control, at their normal maximum operating concentration. Active ingredients tested include: glutaraldehyde, didecyl-dimethyl-ammonium chloride, isolthazoles, dibromonitrilopropionamide, ionic and non-ionic surfactants. The test gives a false positive if used with polymeric biguanide or THPS.
Operator time	5 - 10 minutes
Development time	25 minutes
Sample point	Supplied with both $\frac{1}{2}$ " and $\frac{3}{4}$ " female threaded fittings, additional
	sample points are available (part no. 100171)



hydrosense®

Springfield, Linlithgow EH49 7NX, United Kingdom. T: +44 (0)1506 840888. F: +44 (0)1506 840777.

www.hydrosense.biz

TECH0019v1.0