

Get results in the field. In minutes.

# KATHON™ Test Kit

*An innovative microbicide test kit designed for use with KATHON preservatives from Dow Microbial Control.*

## Easy to use

- Complete the test in 6 steps, 10 minutes
- Perform test in the field or in the lab

## Reliable

- Determine KATHON concentrations via visual colorimetric methods
- Get results with no interference from cooling water additives

## Sensitive

- Measure within 0.75 ppm of the active ingredient
- Test with a range of up to 7.5 ppm for the active ingredient (without dilution)

## Optimize treatment

- Confirm dose amounts and monitor residual concentrations



**DISTRIBUTED IN THE U.S.A. BY**

[www.TheWaterTreatmentStore.com](http://www.TheWaterTreatmentStore.com)